



URINE DRUG SCREEN



One Step Multi-Drug Screen Test Cup (Urine)

All-in-one urine based multi-drug test solutions that can be utilized in practically any drug screening scenario. Our drug test cups can offer various customizable panels with an optional 6 specimen validity tests for urine adulteration.

ROUND CUP (V)

Max. 16 Strips



FLAT CUP (F)

Max. 14 Strips



FEATURES

- Cup V: Clicker audibly confirms when cap is closed securely
- Cup F: Easily and seamlessly accommodates various readers
- Cap ensures secure sealing during transport to laboratories for confirmatory testing
- Health Canada licensed for Class II and Class III (POC + OTC use)
- Class II license includes 81 test items
- Class III license covers 27 test items, including the Canadian clients' most demanded items: 6MAM, EDDP, ETG, FEN (FYL), HMO, KET, MPD, TRA, etc.
- Flexible configuration settings
- French IFU available upon request

SPECIFICATIONS

- Minimum specimen volume: 30 mL
- Qualitative results in 5 minutes
- POC testing or home testing
- Internal, procedural control included
- Specimen validity tests provide up to 6 urine specimen adulteration tests (Oxidants, Specific Gravity, pH, Nitrates, Glutaraldehyde, and Creatine)
- 24-Month shelf life
- Individually packaged cups: 25 tests/kit, 4 kits/shipping carton
- Custom configurations, private-labels, and other cup models available

Healgen Scientific, LLC 3818 Fuqua St. Houston, TX 77047, USA www.healgen.com

Healstone Biotech Inc., Unit #650-655 West Kent Avenue North, Vancouver, B.C., Canada, V6P 6T7

1-800-919-3622/1-604-800-2434/ info.healstone@healgen.com www.healstone.com




TEST MENU

One Step Multi-Drug Screen Test Cup (Urine)



Test	Calibrator	Cut-off (ng/mL)	Class II Multi-Panel Cup	Class III Multi-Panel Cup
Acetaminophen (ACE)	Acetaminophen	5,000	Y	-
Amphetamine (AMP)	D-Amphetamine	250/300	Y	-
		500/1000	Y	Y
Barbiturates (BAR)	Secobarbital	200	Y	-
		300	Y	Y
Buprenorphine (BUP)	Buprenorphine	5	Y	-
		10	Y	Y
Benzodiazepines (BZO)	Oxazepam	150/200	Y	-
		300	Y	Y
Cocaine (COC)	Benzoylecgonine	100	Y	-
		150/300	Y	Y
Cotinine (COT)	Cotinine	200	Y	(Coming soon)
Clonazepam (CLO)	Clonazepam	500	Y	-
MDMA (Ecstasy)	D, L-3,4-Methylenedioxy methamphetamine (MDMA)	150/250/300/2000	Y	-
		500	Y	Y
Ethyl Glucuronide (ETG)	Ethyl Glucuronide	300	Y	-
		500	Y	(Coming soon)
Fentanyl (FEN)	Fentanyl	50/200/300	Y	-
		100	Y	Y
	Nor-fentanyl	10	Y	(Coming soon)
		20	Y	-

TEST PROCEDURE






TEST MENU

One Step Multi-Drug Screen Test Cup (Urine)



Test	Calibrator	Cut-off (ng/mL)	Class II Multi-Panel Cup	Class III Multi-Panel Cup
Gabapentin (GAB)	Gabapentin	1000/2000	Y	-
Hydrocodone (HCD)	Hydrocodone	10/300	Y	-
Hydromorphone (HMO)	Hydromorphone	300	Y	(Coming soon)
K2 Synthetic Cannabinoid	JWH-073/JWH-018	25/50	Y	-
K3 (AB-Pinaca)	AB-Pinaca	10	Y	-
Ketamine (KET)	Ketamine	100	Y	(Coming soon)
		500/1000	Y	-
Kratom (KRA)	Mitragynine	250/300	Y	-
Lysergic acid diethylamide (LSD)	D-lysergic acid diethylamide	10/20	Y	-
Marijuana (THC)	11-nor- Δ^9 -THC-9 COOH	18/20/25/40	Y	-
		50	Y	Y
EDDP (Methadone Metabolites)	2-Ethylidene-1,5-dimethyl-3,3-dipheylpyrrolidine (EDDP)	100	Y	Y
		300	Y	-
Methadone (MTD)	Methadone	50/200	Y	-
		300	Y	Y
Methamphetamine (MET)	D-Methamphetamine	250/300	Y	-
		500/1000	Y	Y
Methaqualone (MQL)	Methaqualone	300	Y	-
Methcathinone (MTC)	Methcathinone	300	Y	-

TEST PROCEDURE



Healgen Scientific, LLC 3818 Fuqua St. Houston, TX 77047, USA www.healgen.com

Healstone Biotech Inc., Unit #650-655 West Kent Avenue North, Vancouver, B.C., Canada, V6P 6T7

1-800-919-3622/1-604-800-2434/ info.healstone@healgen.com www.healstone.com



TEST MENU

One Step Multi-Drug Screen Test Cup (Urine)



Test	Calibrator	Cut-off (ng/mL)	Class II Multi-Panel Cup	Class III Multi-Panel Cup
Methylenedioxypropylamphetamine (MDPV)	3,4-Methylenedioxypropylamphetamine	1000	Y	-
Methylphenidate (MPD)	Methylphenidate	300	Y	(Coming soon)
6-Monoacetylmorphine (6-MAM)	6-Monoacetylmorphine	10	Y	(Coming soon)
		20	Y	-
Morphine (MOP 300)	Morphine	300	Y	Y
Morphine (MOP 150)		150	Y	-
Morphine (OPI, MOP 2000)		2000	Y	Y
Oxycodone (OXY)	Oxycodone	100	Y	Y
Phencyclidine (PCP)	Phencyclidine	25	Y	Y
Pregabalin (PGB)	Pregabalin	500/1000	Y	-
Propoxyphene (PPX)	Propoxyphene	300	Y	-
UR-144	Synthetic Cannabinoid	50	Y	-
Carisoprodol (SOMA)	Carisoprodol	1000	Y	-
Tramadol (TRA)	Tramadol	100/300	Y	-
		200	Y	(Coming soon)
Tricyclic Antidepressants (TCA)	Nortriptyline	300	Y	-
		1000	Y	Y
Alcohol (ALC)	Ethanol	>0.04%B.A.C	Y	-
Oxidants (OX)	-	0.30%	Y	(Coming soon)
Specific Gravity (S.G.)	-	0.21%	Y	(Coming soon)
pH	-	0.06%	Y	(Coming soon)
Nitrite (NIT)	-	0.06%	Y	(Coming soon)
Glutaraldehyde (GLU)	-	0.02%	Y	(Coming soon)
Creatinine (CRE)	-	0.03%	Y	(Coming soon)