

HEALSTONE BIOTECH INC.

Product Data: Fentanyl (FEN) BSA Antigen (Low Hapten)

Product Type	Antigen
Product Name	FEN-BSA (L)
Catalog No.	FEN-HS-Test
Description	BSA conjugated Fentanyl antigen with low hapten
Form	Liquid
Supplied In	Phosphate Buffered Saline pH 7.2
Preservative	0.02% Sodium Azide
Applications	Lateral Flow, ELISA
Matched Pair	Pairs with most FEN AB on the market
Storage	Short term (≤ 2 weeks): 2-8°C. Long term: -20°C. Avoid repeated freeze-thaw cycles.
Product Category	Drugs of Abuse
Note	This product is intended for research and manufacturing uses only. It is not for human, therapeutic, or in-vivo use. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. Made in Canada.

