

HEALSTONE BIOTECH INC.

Product Data:

Synthetic Cannabinoids (K3, AB-Pinaca) Monoclonal Antibody

| | |
|------------------|--|
| Product Type | Antibody |
| Product Name | K3MAB |
| Catalog No. | K3MAB -HS-001 |
| Description | Chromatography purified anti-AB-Pinaca monoclonal antibody |
| Immunogen | AB-Pinaca-BSA |
| Host/Source | Mouse |
| Specificity | AB-Pinaca; AB-Pinaca 4-Hydroxypentyl metabolite; AB-Pinaca 5-Hydroxypentyl metabolite; AB-Pinaca 5-Pentanoic acid metabolite |
| Form | Liquid |
| Supplied In | Phosphate Buffered Saline pH 7.2 |
| Preservative | 0.05% Sodium Azide |
| Matched Pair | K3-HS-001 |
| Applications | Lateral Flow, ELISA |
| Storage | Short term: 2-8°C. Long term: -20°C. Avoid repeated freeze-thaw cycles. |
| Product Category | Drugs of Abuse |
| Note | This product is intended for research and manufacturing uses only. It is not for human, therapeutic, or in-vivo use. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. Made in Canada. |

