

HEALSTONE BIOTECH INC.

Product Data: Lysergic acid diethylamide (LSD) Monoclonal Antibody

Product Type	Antibody
Product Name	LSDMAB
Catalog No.	LSDMAB-HS-002
Description	Chromatography purified anti-Lysergic Acid Diethylamide monoclonal antibody
Immunogen	Lysergic Acid Diethylamide-KLH
Host/Source	Mouse
Specificity	Lysergic Acid Diethylamide (LSD)
Form	Liquid
Supplied In	Phosphate Buffered Saline pH 7.2
Preservative	0.05% Sodium Azide
Matched Pair	LSD-HS-001
Applications	Lateral Flow, ELISA
Storage	Short term: 2-8°C. Long term: -20°C. Avoid repeated freeze-thaw cycles.
Product Category	Drugs of Abuse
Note	This product is intended for research and manufacturing uses only. It is not for human, therapeutic, or in-vivo use. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. Made in Canada.

