

# HEALSTONE BIOTECH INC.

---

## Product Data: Methadone (MTD) BSA Antigen

Product Type	Antigen
Product Name	MTD-BSA
Catalog No.	MTD-HS-001
Description	BSA conjugated Methadone antigen
Form	Liquid
Supplied In	Phosphate Buffered Saline pH 7.2
Preservative	0.02% Sodium Azide
Applications	Lateral Flow, ELISA
Storage	Short term ( $\leq 2$ weeks): 2-8°C. Long term: -20°C. Avoid repeated freeze-thaw cycles.
Product Category	Drugs of Abuse
Note	This product is intended for research and manufacturing uses only. It is not for human, therapeutic, or in-vivo use. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. Made in Canada.

