

# HEALSTONE BIOTECH INC.

---

## Product Data: Phencyclidine (PCP) BSA Antigen

|                  |  |
|------------------|--|
| Product Type     | Antigen  |
| Product Name     | PCP-BSA  |
| Catalog No.      | PCP-HS-001   |
| Description      | BSA conjugated Phencyclidine antigen   |
| Form             | Liquid   |
| Supplied In      | Phosphate Buffered Saline pH 7.2   |
| Preservative     | 0.02% Sodium Azide   |
| Applications     | Lateral Flow, ELISA  |
| Storage          | Short term ( $\leq 2$ weeks): 2-8°C. Long term: -20°C.<br>Avoid repeated freeze-thaw cycles.   |
| Product Category | Drugs of Abuse   |
| Note             | This product is intended for research and manufacturing uses only. It is not for human, therapeutic, or in-vivo use. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. Made in Canada. |

