

HEALSTONE BIOTECH INC.

Product Data: 25-Hydroxy Vitamin D3 (VD3) BSA Antigen

| | |
|------------------|--|
| Product Type | Antigen |
| Product Name | VD3-BSA |
| Catalog No. | VD3-HS-001 |
| Description | BSA conjugated 25-Hydroxy Vitamin D3 antigen |
| Form | Liquid |
| Supplied In | Phosphate Buffered Saline pH 7.2 |
| Preservative | 0.02% Sodium Azide |
| Applications | Lateral Flow, ELISA |
| Storage | Short term (≤ 2 weeks): 2-8°C. Long term: -20°C. Avoid repeated freeze-thaw cycles. |
| Product Category | Drugs of Abuse |
| Note | This product is intended for research and manufacturing uses only. It is not for human, therapeutic, or in-vivo use. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. Made in Canada. |

