

# HEALSTONE BIOTECH INC.

---

## Product Data:

### Goat anti-Human Serum Albumin (hALB) Polyclonal Antibody

Product Type	Antibody
Product Name	hALBPAB
Catalog No.	hALBPAB -HS-001
Description	Chromatography purified anti-human serum albumin polyclonal antibody
Immunogen	Human albumin
Host/Source	Goat
Specificity	Human Albumin, does not cross-react with mouse and bovine serum albumin.
Form	Liquid
Supplied In	Phosphate Buffered Saline pH 7.2
Preservative	0.02% Sodium Azide
Applications	Lateral Flow, ELISA
Storage	Short term: 2-8°C. Long term: -20°C. Avoid repeated freeze-thaw cycles.
Product Category	Drugs of Abuse
Note	This product is intended for research and manufacturing uses only. It is not for human, therapeutic, or in-vivo use. It is not a diagnostic device. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations. Made in Canada.

