

## Rapid Test Cassettes

Suitable for reagent strip loading and test window information display for immunochromatographic reactions

- \* Single and multiple test windows available
- \* Multiple models and sizes available
- \* Custom printing available
- \* High-quality PS material ensures hardness while retaining toughness.



HSRTC-1



HSRTC-2



HSRTC-3



HSRTC-5



HSRTC-6



HSRTC-7



HSRTC-8



HSRTC-9



HSRTC-10



HSRTC-11



HSRTC-12




HSRTC-13




HSRTC-14




HSRTC-15



HSRTC-16



HSRTC-17



HSRTC-18



HSRTC-19



Cat. No.	P/N	Material	Specification	Packaging
HSRTC-1	Cassette	PS	3.95 mm * 60 mm	
HSRTC-2	Cassette	PS	4.00 mm * 60 mm	
HSRTC-3	Cassette	PS	3.00 mm * 60 mm	
HSRTC-4	Cassette	PS	3.95 mm * 60 mm	200 pcs/bag, 20 bags/cs
HSRTC-5	Cassette	PS	3.95 mm * 60 mm	
HSRTC-6	Cassette	PS	3.95 mm * 60 mm	
HSRTC-7	Cassette	PS	3.50 mm * 60 mm	
HSRTC-8	Cassette	PS	4.05 mm * 60 mm	150 pcs/bag, 20 bags/cs
HSRTC-9	Cassette	PS	3.00 mm * 60 mm	400 pcs/bag, 20 bags/cs
HSRTC-10	Cassette	PS	4.25 mm * 60 mm	150 pcs/bag, 20 bags/cs
HSRTC-11	Cassette	PS	4.25 mm * 60 mm	150 pcs/bag, 20 bags/cs
HSRTC-12	Cassette	PS	4.10 mm * 60 mm	100 pcs/bag, 20 bags/cs
HSRTC-13	Cassette	PS	4.00 mm * 70 mm	180 pcs/bag, 20 bags/cs
HSRTC-14	Cassette	PS	3.93 mm * 74 mm	50 pcs/bag, 20 bags/cs
HSRTC-15	Cassette	PS	4.00 mm * 60 mm	125 pcs/bag, 20 bags/cs
HSRTC-16	Cassette	PS	4.00 mm * 60 mm	125 pcs/bag, 20 bags/cs