



# Avian Influenza Virus (AIV) (H5) Antigen Rapid Tests



The Healgen® Avian Influenza Virus (AIV) (H5) Antigen Rapid Test is a rapid, lateral flow immunochromatographic assay intended for the qualitative, differential detection of H5 subtypes of AIV virus hemagglutinin protein in samples such as nasal swab, deep throat/oropharyngeal swab, anal/cloacal swab, feces and milk specimens from animals suspected of AIV Infections.

(For Research Use Only)

## FEATURES

- Rapid testing for AIV H5 subtypes antigen within 15 minutes
- Specimen include nasal swab, throat/oropharyngeal swab, anal/cloacal swab, fecal and milk
- Supports quick and informed treatment decisions in veterinary practice
- All necessary reagents provided & no equipment needed.
- High sensitivity and no cross-reactivity with other Flu A subtypes

## SPECIFICATIONS

- 10 tests per kit
- Storage: 2-30°C/36-86°F
- Shelf life: 24 months






# Avian Influenza Virus (AIV) (H5) Antigen Rapid Tests



## Test Procedures (Reference only, please strictly follow IFU):



## Ordering Information

Product Description	Specimen	Catalog No.	Kit Size
Avian Influenza Virus (AIV) (H5) Antigen Rapid Test	nasal swab, throat/ oropharyngeal swab, anal/cloacal swab, fecal and milk	GCFLU(AH5)-1702a	10 Tests/Kit