



# Influenza A Virus (Flu A) Antigen Rapid Tests



The Healgen® Influenza A Virus (Flu A) Antigen Rapid Test is a rapid, lateral flow immunochromatographic assay intended for the qualitative, rapid differential detection of influenza A virus nucleocapsid protein (NP) in samples including nasal swab, throat/oropharyngeal swab, anal/cloacal swab, fecal and milk specimen from animals suspected of IAV Infections.

(For Research Use Only)

## FEATURES

- Rapid testing for Flu A antigen within 15 minutes
- Specimen type include nasal swab, throat/oropharyngeal swab, anal/cloacal swab, fecal and milk
- Quick procedure supports informed treatment decisions in veterinary practice
- All necessary reagents provided & no equipment needed
- High sensitivity and no cross-reactivity with leading respiratory microbes
- Inclusivity spanning various strains of Influenza A
- Developed using regionally and genetically diverse sources of viral samples



## SPECIFICATIONS


- 10 tests per kit
- Storage: 2-30°C/36-86°F
- Shelf life: 24 months



# Influenza A Virus (Flu A) Antigen Rapid Tests



## Test Procedures (Reference only, please strictly follow IFU):



## Ordering Information

Product Description	Specimen	Catalog No.	Kit Size
Influenza A Virus (Flu A) Antigen Rapid Test	nasal swab, throat/ oropharyngeal swab, anal/ cloacal swab, fecal and milk	GCFLU(A)-1702a	10 Tests/Kit