



**Healstone
Biotech**

Products in the IVD World

Raw Materials

2025



ABOUT US

Founded in 2016 by Canadian researchers in Vancouver, Healstone Biotech Inc. is at the forefront of biotechnology, dedicated to the advancement of diagnostic solutions. Our commitment to innovation and excellence has helped forged strong partnerships with international clientele across the United States, Canada, Europe, and China. At Healstone Biotech, we specialize in:

- Development of advanced immunological reagents
- Research, design, and development of IVD rapid test formats
- Distribution of a wide range of IVD rapid tests
- Supplying premium laboratory materials used for IVD development and other academic or industry applications

At Healstone Biotech Inc., we are proud to offer a wide range of exceptional raw materials for lateral flow assay development applications, including synthesized and native antigens, primary and secondary antibodies, and control materials. Our products are essential components for the development of rapid diagnostic tests, tailored to fulfill needs in various categories including toxicology, infectious diseases, cardiovascular, animal diseases, and others. Choose us as your trusted partner to support your diagnostic solutions and join a community of forward-thinking scientists dedicated to improving global healthcare.



TABLE OF CONTENTS

- 01 Toxicology
- 02 Cardiovascular
- 03 Tumor Markers
- 04 Infectious Diseases
- 05 Inflammation
- 06 Glucose Metabolism
- 07 Endocrine
- 08 Human Serum
- 09 Secondary Antibodies
- 10 Animal Diseases
- 11 Blocking Reagents
- 12 Miscellaneous



01. Toxicology

* mAb = Monoclonal antibody

pAb = Polyclonal antibody

Ag = Antigen

Analyte	Product *	Cat. #	Source	Applications in LFA
25-Hydroxy Vitamin D3 (VD3)	VD3-BSA Ag	VD3-HS-001	Synthetic	Coating
7-Aminoclonazepam (7-ACL)	7-ACL mAb	HS-A4125	Mouse	Labeling
	7-ACL-BSA Ag	HS-A4126	Synthetic	Coating
Alprazolam (ALP)	ALP mAb	HS-A4133	Mouse	Labeling
	ALP-BSA Ag	HS-A4134	Synthetic	Coating
Amphetamine (AMP)	AMP-BSA Ag	AMP-HS-001	Synthetic	Coating
	AMP mAb	HS-A4013	Mouse	Labeling
	AMP-BSA Ag	HS-A4088	Synthetic	Coating
Barbiturates (BAR)	BAR-BSA Ag	BAR-HS-001	Synthetic	Coating
	BAR mAb	HS-A4015	Mouse	Labeling
	BAR-BSA Ag	HS-A4121	Synthetic	Coating
Benzodiazepine (BZO)	BZO mAb	HS-A4019	Mouse	Labeling
	BZO-BSA Ag	HS-A4092	Synthetic	Coating
Buprenorphine (BUP)	BUP-BSA Ag	BUP-HS-001	Synthetic	Coating
	BUP mAb	HS-A4017	Mouse	Labeling
	BUP-BSA Ag	HS-A4096	Synthetic	Coating
Caffeine (CAF)	CAF-BSA Ag	CAF-HS-001	Synthetic	Coating
	CAF mAb	HS-A4153	Mouse	Labeling
	CAF-BSA Ag	HS-A4154	Synthetic	Coating
Carfentanil (CFYL)	CFYL mAb	HS-A4148	Mouse	Labeling
	CFYL-BSA Ag	HS-A4149	Synthetic	Coating
Carisoprodol SOMA (SOM)	SOMA pAb	SOMPAB-HS-001	Rabbit	Labeling
	SOMA-BSA Ag	SOM-HS-001	Synthetic	Coating
	SOMA mAb	HS-A4100	Mouse	Labeling
	SOMA-BSA Ag	HS-A4101	Synthetic	Coating
Cocaine (COC)	COC-BSA Ag	COC-HS-001	Synthetic	Coating
	COC mAb	HS-A4021	Mouse	Labeling
	COC-BSA Ag	HS-A4022	Synthetic	Coating
Codeine (COD)	COD mAb	HS-A4123	Mouse	Labeling
	COD-BSA Ag	HS-A4124	Synthetic	Coating
Cotinine, Nicotine Metabolite (COT)	COT-BSA Ag	COT-HS-001	Synthetic	Coating
	COT mAb	HS-A4179	Mouse	Labeling
	COT-BSA Ag	HS-A4118	Synthetic	Coating
Ethyl Glucuronide, Ethanol Metabolite (ETG)	ETG-BSA Ag	ETG-HS-001	Synthetic	Coating
	ETG mAb	HS-A4043	Goat	Labeling
	ETG-BSA Ag	HS-A4087	Synthetic	Coating



01. Toxicology

Analyte	Product *	Cat. #	Source	Applications in LFA
Fentanyl (FEN)	FEN-BSA Ag (High Molar Ratio)	FEN-HS-001	Synthetic	Coating
	FEN-BSA Ag (Low Molar Ratio)	FEN-HS-Test	Synthetic	Coating
	FEN mAb	HS-A4029	Mouse	Labeling
	FEN-BSA Ag	HS-A4098	Synthetic	Coating
Flaka (α -PVP)	α -PVP mAb	HS-A4169	Mouse	Labeling
	α -PVP-BSA Ag	HS-A4170	Synthetic	Coating
Fluoroamine Ketone (F-KET)	F-KET mAb	HS-A4150	Mouse	Labeling
		HS-A4173	Mouse	Labeling
	F-KET-BSA Ag	HS-A4068	Synthetic	Coating
Fluoxetine (FLX)	FLX mAb	HS-A4144	Mouse	Labeling
	FLX-BSA Ag	HS-A4145	Synthetic	Coating
Gabapentin (GAB)	GAB mAb	HS-A4108	Mouse	Labeling
	GAB-BSA Ag	HS-A4109	Synthetic	Coating
Hydrocodone (HCD)	HCD mAb	HS-A4131	Mouse	Labeling
	HCD-BSA Ag	HS-A4132	Synthetic	Coating
Ketamine (KET)	KET mAb	KETMAB-HS-001	Mouse	Labeling
	KET-BSA Ag	KET-HS-001	Synthetic	Coating
	KET mAb	HS-A4041	Mouse	Labeling
	KET-BSA Ag	HS-A4042	Synthetic	Coating
Kratom, Mitragynine (KRA)	KRA mAb	KRAMAB-HS-002	Mouse	Labeling
	KRA-BSA Ag	KRA-HS-003H	Synthetic	Coating
	KRA-HRP Ag	KRA-HS-004	Synthetic	Coating
	KRA mAb	HS-A4110	Mouse	Labeling
	KRA-BSA Ag	HS-A4111	Synthetic	Coating
Lysergic Acid Diethylamide (LSD)	LSD mAb	LSDMAB-HS-002	Mouse	Labeling
	LSD-BSA Ag	LSD-HS-001	Synthetic	Coating
4-methylmethcathinone, Mephedrone (MCAT, CAT)	MCAT-BSA Ag	CAT-HS-001	Synthetic	Coating
3,4-Methylenedioxy Methamphetamine (MDMA)	MDMA-BSA Ag	MDMA-HS-001	Synthetic	Coating
	MDMA mAb	HS-A4001	Mouse	Labeling
	MDMA-BSA Ag	HS-A4095	Synthetic	Coating



01. Toxicology

Analyte	Product *	Cat. #	Source	Applications in LFA
Mescaline (MES)	MES-BSA Ag	MES-HS-001	Synthetic	Coating
	MES mAb	HS-A4183	Mouse	Labeling
	MES-BSA Ag	HS-A4184	Synthetic	Coating
Meperidine, Demerol (MPRD)	MPRD mAb	HS-A4164	Mouse	Labeling
	MPRD-BSA Ag	HS-A4165	Synthetic	Coating
Methadone (MTD)	MTD-BSA Ag	MTD-HS-001	Synthetic	Coating
	MTD mAb	HS-A4005	Mouse	Labeling
	MTD-BSA Ag	HS-A4006	Synthetic	Coating
Methadone metabolites (EDDP)	EDDP-BSA Ag	EDDP-HS-001	Synthetic	Coating
	EDDP mAb	HS-A4025	Mouse	Labeling
	EDDP-BSA Ag	HS-A4026	Synthetic	Coating
Methamphetamine (MET)	MET-BSA Ag	MET-HS-001	Synthetic	Coating
	MET mAb	HS-A4007	Mouse	Labeling
	MET-BSA Ag	HS-A4116	Synthetic	Coating
Methaqualone (MQL)	MQL-BSA Ag	MQL-HS-001	Synthetic	Coating
	MQL mAb	HS-A4151	Mouse	Labeling
	MQL-BSA Ag	HS-A4152	Synthetic	Coating
Methylenedioxypropylvalerone (MDPV)	MDPV mAb	HS-A4185	Mouse	Labeling
	MDPV-BSA Ag	HS-A4186	Synthetic	Coating
Methylphenidate (MPD)	MPD-BSA Ag	MPD-HS-001	Synthetic	Coating
Morphine (MOP)	MOP mAb	MOPMAB-HS-001	Mouse	Labeling
	MOP-BSA Ag	MOP-HS-001	Synthetic	Coating
	MOP mAb	HS-A4011	Mouse	Labeling
	MOP-BSA Ag	HS-A4012	Synthetic	Coating
Opioids (OPI)	OPI-BSA Ag	HS-A4053	Synthetic	Labeling
	OPI mAb	HS-A4052	Mouse	Coating
Oxycodone (OXY)	OXY-BSA Ag (High molar ratio)	OXY-HS-001	Synthetic	Coating
	OXY-BSA Ag (Low molar ratio)	OXY-HS-002	Synthetic	Coating
	OXY mAb	HS-A4003	Mouse	Labeling
	OXY-BSA Ag	HS-A4004	Synthetic	Coating
Phenyl Cyclo Piperidine (PCP)	PCP-BSA Ag	PCP-HS-001	Synthetic	Coating
	PCP mAb	HS-A4039	Mouse	Labeling
	PCP-BSA Ag	HS-A4040	Synthetic	Coating




01. Toxicology

Analyte	Product *	Cat. #	Source	Applications in LFA
Phenytoin (PYT)	PYT-BSA Ag	PYT-HS-001	Synthetic	Coating
Pregabalin (PGB)	PGB mAb	HS-A4106	Mouse	Labeling
	PGB-BSA Ag	HS-A4107	Synthetic	Coating
Propoxyphene (PPX)	PPX mAb	HS-A4033	Mouse	Labeling
	PPX-BSA Ag	HS-A4120	Synthetic	Coating
Quetiapine (QTP)	QTP mAb	HS-A4138	Mouse	Labeling
	QTP-BSA Ag	HS-A4139	Synthetic	Coating
Risperidone (RPD)	RPD mAb	HS-A4140	Mouse	Labeling
	RPD-BSA Ag	HS-A4141	Synthetic	Coating
Synthetic Cannabinoids K2	K2 mAb	HS-A4037	Mouse	Labeling
	K2-BSA Ag	HS-A4097	Synthetic	Coating
Synthetic Cannabinoids K3 (AB-PINACA)	K3 mAb	K3MAB-HS-001	Mouse	Labeling
	K3-BSA Ag	K3-HS-001	Synthetic	Coating
	K3 mAb	HS-A4102	Mouse	Labeling
	K3-BSA Ag	HS-A4103	Synthetic	Coating
Synthetic Cannabinoids K4 (UR-144)	K4 mAb	HS-A4104	Mouse	Labeling
	K4-BSA Ag	HS-A4105	Synthetic	Coating
Tetrahydrocannabinol (THC)	THC mAb	HS-A4027	Mouse	Labeling
	THC-BSA Ag	HS-A4114	Synthetic	Coating
Tramadol (TRA)	TRA mAb	TRAMAB-HS-001	Mouse	Labeling
	TRA-BSA Ag	TRA-HS-001	Synthetic	Coating
	TRA mAb	HS-A4035	Mouse	Labeling
	TRA-BSA Ag	HS-A4112	Synthetic	Coating
Trazodone (TZD)	TZD mAb	HS-A4135	Mouse	Labeling
	TZD-BSA Ag	HS-A4136	Synthetic	Coating
Tricyclic Antidepressant (TCA)	TCA-BSA Ag	TCA-HS-001	Synthetic	Coating
	TCA mAb	HS-A4009	Mouse	Labeling
	TCA-BSA Ag	HS-A4189	Synthetic	Coating
Tropicamide (TRO)	TRO mAb	HS-A4142	Mouse	Labeling
	TRO-BSA Ag	HS-A4143	Synthetic	Coating
Xylazine (XYL)	XYL mAb	HS-A4175	Mouse	Labeling
	XYL-BSA Ag	HS-A4188	Synthetic	Coating
Zolpidem (ZOL)	ZOL mAb	HS-A4129	Mouse	Labeling
	ZOL-BSA Ag	HS-A4130	Synthetic	Coating
Zopiclone (ZOP)	ZOP mAb	HS-A4127	Mouse	Labeling
	ZOP-BSA Ag	HS-A4128	Synthetic	Coating



02. Cardiovascular

Analyte	Product *	Cat. #	Source	Applications in LFA
Cardiac Troponin I (cTnI)	cTnI mAb	HS-A3001	Mouse	Coating
		HS-A3002	Mouse	Labeling
		HS-A3024	Mouse	Labeling
		HS-A3025	Mouse	Coating
	hs-cTnI mAb	HS-A3049	Mouse	Labeling
		HS-A3050	Mouse	Coating
		HS-A3051	Mouse	Labeling
Cardiac Troponin T (cTnT)	cTnT mAb	HS-A3016	Mouse	Coating
		HS-A3017	Mouse	Labeling
		HS-A3026	Mouse	Labeling
		HS-A3027	Mouse	Coating
Copeptin (CPP)	Recombinant CPP Ag	HS-M1026	E. coli	Control Material
	CPP mAb	HS-M1020	Mouse	Labeling
		HS-M1021	Mouse	Coating
Creatine Kinase Isoenzyme-MB (CK-MB)	CK-MB mAb	HS-A3007	Mouse	Coating
		HS-A3008	Mouse	Labeling
		HS-A3028	Mouse	Labeling
		HS-A3029	Mouse	Coating
		HS-A3053	Mouse	—
		HS-A3054	Mouse	—
D-Dimer (D-D)	D-D mAb	HS-A3020	Mouse	Coating
		HS-A3021	Mouse	Labeling
		HS-A3032	Mouse	Labeling
		HS-A3033	Mouse	Coating
		HS-A3037	Mouse	Labeling
		HS-A3038	Mouse	Coating





02. Cardiovascular

Analyte	Product *	Cat. #	Source	Applications in LFA
Growth Stimulation Expressed Gene 2 (ST2)	Recombinant ST2 Ag	HS-A3046	293 cells	Control Material
	ST2 mAb	HS-A3047	Mouse	Labeling
		HS-A3048	Mouse	Coating
Myoglobin (Myo)	Myo mAb	HS-A3004	Mouse	Coating
		HS-A3005	Mouse	Labeling
		HS-A3030	Mouse	Labeling
		HS-A3031	Mouse	Coating, Labeling
		HS-A3040	Mouse	Coating
N-terminal Brain Natriuretic Peptide Precursor (NT proBNP)	NT proBNP mAb	HS-A3022	Mouse	Coating
		HS-A3023	Mouse	Labeling
		HS-A3034	Mouse	Labeling
		HS-A3035	Mouse	Coating
		HS-A3039	Mouse	Labeling
Oxidized Low Density Lipoprotein (OxLDL)	OxLDL Ag	HS-M1027	Native	Control Material
	OxLDL mAb	HS-M1022	Mouse	Labeling
		HS-M1023	Mouse	Coating
S100 β Protein (S100 β)	Recombinant S100 β Ag	HS-M1028	E. coli	Control Material
		HS-M1036	E. coli	Control Material
	S100 β mAb	HS-M1024	Mouse	Labeling
		HS-M1025	Mouse	Coating



Fertility & Women's Health



Drugs



Infectious Diseases



Oncology



Cardiovascular



Autoimmunity



Metabolism



03. Tumor Markers

Analyte	Product *	Cat. #	Source	Applications in LFA
Alpha-Fetoprotein (AFP)	AFP mAb	HS-A2007	Mouse	Labeling
		HS-A2008	Mouse	Coating
		HS-A2044	Mouse	Labeling
		HS-A2045	Mouse	Coating
		HS-A2054	Mouse	Coating
Carbohydrate Antigen 125 (CA125)	CA125 mAb	HS-A2022	Mouse	Labeling
		HS-A2023	Mouse	Coating
		HS-A2073	Mouse	Labeling
		HS-A2074	Mouse	Coating
Carbohydrate Antigen 153 (CA153)	CA153 mAb	HS-A2047	Mouse	Labeling
		HS-A2048	Mouse	Coating
Carbohydrate Antigen 199 (CA199)	CA199 mAb	HS-A2034	Mouse	Coating
		HS-A2035	Mouse	Labeling
Carbohydrate Antigen 242 (CA242)	CA242 mAb	HS-A2086	Mouse	Labeling
		HS-A2087	Mouse	Coating
Carcinoembryonic Antigen (CEA)	CEA mAb	HS-A2005	Mouse	Labeling
		HS-A2006	Mouse	Coating
Cytokeratin (Cy 211)	Cy 211 mAb	HS-A2029	Mouse	Coating
		HS-A2046	Mouse	Labeling
		HS-A2071	Mouse	Labeling
Fecal Occult Blood (FOB)	FOB mAb	HS-A2001	Mouse	Labeling
		HS-A2002	Mouse	Coating
		HS-A2025	Mouse	Labeling
		HS-A2026	Mouse	Coating
Ferritin (FER)	FER mAb	HS-A2040	Mouse	Labeling
		HS-A2041	Mouse	Coating
Gastrin 17 (G17)	Recombinant G17 Ag	HS-M1029	E. coli	Control Material
	G17 mAb	HS-M1018	Mouse	Labeling
		HS-M1019	Mouse	Coating
		HS-M1034	Mouse	Coating
		HS-M1035	Mouse	Labeling
Human Epididymis Protein 4 (HE4)	HE4 mAb	HS-A2036	Mouse	Labeling
		HS-A2037	Mouse	Coating
		HS-A2060	Mouse	Labeling
		HS-A2085	Mouse	Labeling



03. Tumor Markers

Analyte	Product *	Cat. #	Source	Applications in LFA
Neuron-Specific Enolase (NSE)	NSE mAb	HS-A2038	Mouse	Labeling
		HS-A2072	Mouse	Coating
Pepsinogen (PG)	PG I mAb	HS-A2049	Mouse	Labeling
		HS-A2050	Mouse	Coating
		HS-A2061	Mouse	Labeling
		HS-A2062	Mouse	Coating
		HS-A2063	Mouse	Labeling
		HS-A2064	Mouse	Coating
	PG II mAb	HS-A2051	Mouse	Labeling
		HS-A2052	Mouse	Coating
		HS-A2065	Mouse	Labeling
		HS-A2066	Mouse	Coating
		HS-A2067	Mouse	Labeling
		HS-A2068	Mouse	Coating
Pro-Gastrin-Releasing Peptide (ProGRP)	ProGRP mAb	HS-A2027	Mouse	Coating
		HS-A2028	Mouse	Labeling
Prostate Specific Antigen (PSA)	PSA mAb	HS-A2003	Mouse	Labeling
		HS-A2004	Mouse	Coating
Squamous Cell Carcinoma (SCCA)	SCCA mAb	HS-A2019	Mouse	Labeling
		HS-A2020	Mouse	Coating
		HS-A2069	Mouse	Labeling
		HS-A2070	Mouse	Coating
Vascular Endothelial Growth Factor (VEGF)	Recombinant VEGF Ag	HS-A2059	E. coli	Control Material
	VEGF mAb	HS-A2055	Mouse	Labeling
		HS-A2056	Mouse	Labeling
		HS-A2057	Mouse	Labeling
	VEGF pAb	HS-A2058	Goat	Coating
HS-A2096		Mouse	Coating	



04. Infectious Diseases

Analyte	Product *	Cat. #	Source	Applications in LFA
Adenovirus (ADV)	ADV mAb	HS-A5045	Mouse	Labeling
		HS-A5046	Mouse	Coating
Chikungunya Virus (CHIK)	Recombinant CHIK Ag	HS-A5125	E. coli	Coating
Chlamydia Pneumoniae (Cpn)	Recombinant Cpn Ag	HS-A5100	E. coli	Coating, Labeling
COVID-19	Recombinant 2019-nCoV NP Ag	HS-A5056	293 cells	Control Material
	Recombinant 2019-nCoV S Ag (Trimer)	HS-A5062	293 cells	Coating, Labeling
	Recombinant 2019-nCoV S1 Ag	HS-A5065	293 cells	Coating, Labeling
	Recombinant hACE2 hFc Ag	HS-A5102	293 cells	Coating, Labeling
	SARS-COV2RBD mAb	HS-A5082	—	Control Material
	2019-nCoV NP mAb	HS-A5089	Mouse	Labeling
		HS-A5090	Mouse	Coating
		HS-A5091	Mouse	Labeling
		HS-A5092	Mouse	Coating
		HS-A5109	Mouse	Coating
		HS-A5134	Rabbit	Labeling
Dengue Virus (DENV)	Recombinant DENV E Ag	HS-A5114	E. coli	Coating
		HS-A5115	E. coli	Coating
Enterovirus EV71	Recombinant EV71-VP1 Ag	HS-M1011	E. coli	Coating, Labeling
Epstein-Barr (EB) virus	Recombinant EB VCA Ag	HS-M1001	E. coli	Coating, Labeling
	Recombinant EB EA Ag	HS-M1002	E. coli	Coating, Labeling
	Recombinant EB NA Ag	HS-M1003	E. coli	Coating, Labeling
	Recombinant EB Rta Ag	HS-M1004	E. coli	Coating, Labeling
Group A Streptococcal Infection DNaseB	Recombinant DNaseB Ag	HS-M1010	E. coli	Coating, Labeling



04. Infectious Diseases

Analyte	Product *	Cat. #	Source	Applications in LFA
Helicobacter pylori (Hp)	Hp Ag	HS-M1007	43504	—
	Recombinant Hp CagA Ag	HS-A5116-2	E. coli	Coating
		HS-A5116-3	E. coli	Coating
	Recombinant Hp VacA Ag	HS-A5120	E. coli	Coating
	Recombinant Hp CagA+VacA Ag	HS-A5121	E. coli	Coating
		HS-A5121-1	E. coli	Coating
		HS-A5121-2	E. coli	Coating
	Hp mAb	HS-A5049	Mouse	Labeling
HS-A5050		Mouse	Coating	
Hepatitis B Core Antigen (HBcAg)	Recombinant HBcAg Ag	HS-A5184	E. coli	Coating, Labeling
Hepatitis B Surface Antigen (HBsAg)	HBsAg Ag	HS-A5003	Native	Coating
	Recombinant HBsAg Ag	HS-A5004	—	Labeling
	HBsAg mAb	HS-A5001	Mouse	Labeling
		HS-A5002	Mouse	Coating
Hepatitis C (HCV)	Recombinant HCV Ag	HS-A5010	E. coli	Coating, Labeling
Hepatitis E Virus (HEV)	Recombinant HEV 1+4 Ag	HS-M1005	E. coli	Coating, Labeling
Human Immunodeficiency Virus (HIV)	Recombinant HIV 1+2 Ag	HS-A5053	E. coli	Labeling
		HS-A5199	E. coli	Coating
	Recombinant HIV gp24 Ag	HS-A5066	E. coli	Control Material
	Recombinant HIV gp36 Ag	HS-A5012	E. coli	Coating
		HS-A5197	E. coli	Labeling
	Recombinant HIV gp41 Ag	HS-A5013	E. coli	Coating
		HS-A5198	E. coli	Labeling
	HIV gp24 mAb	HS-A5015	Mouse	Labeling
		HS-A5016	Mouse	Coating
HS-A5054		Mouse	Labeling	
HS-A5055		Mouse	Coating	
Human Parvovirus B19 (B19)	Recombinant B19 Ag	HS-M1009	E. coli	Coating, Labeling



04. Infectious Diseases

Analyte	Product *	Cat. #	Source	Applications in LFA
Influenza	Recombinant FluA NP Ag	HS-A5147	E. coli	Control Material
	FluA mAb	HS-A5041	Mouse	Labeling
		HS-A5042	Mouse	Coating
		HS-A5180	Mouse	Labeling
		HS-A5181	Mouse	Coating
	Recombinant FluB NP Ag	HS-A5148	E. coli	Control Material
	FluB mAb	HS-A5043	Mouse	Labeling
		HS-A5044	Mouse	Coating
Malaria	Recombinant Mala pf Ag	HS-A5123	E. coli	Coating, Labeling
	Recombinant Mala pv Ag	HS-A5124	E. coli	Coating, Labeling
	Malaria pf mAb	HS-A5032	Mouse	Labeling
		HS-A5033	Mouse	Coating
	Malaria pv mAb	HS-A5030	Mouse	Labeling
		HS-A5031	Mouse	Coating
	Malaria pan mAb	HS-A5034	Mouse	Labeling
		HS-A5035	Mouse	Coating
Mycoplasma Pneumoniae (MP)	MP Ag	HS-A5101	Native	Coating, Labeling
		HS-A5179	Native	Control Material
	Recombinant MP Ag	HS-M1008	E. coli	Coating, Labeling
	MP mAb	HS-A5117	Mouse	Coating
		HS-A5118	Mouse	Labeling
Respiratory Syncytial Virus (RSV)	Recombinant RSV Ag	HS-A5099	E. coli	Control Material
	RSV mAb	HS-A5097	Mouse	Labeling
		HS-A5098	Mouse	Coating
		HS-A5171	Mouse	Labeling
		HS-A5172	Mouse	Coating
Rotavirus (ROTA)	ROTA mAb	HS-A5047	Mouse	Labeling
		HS-A5048	Mouse	Coating
Treponema Pallidum (TP)	Recombinant TP Ag	HS-A5018	E. coli	Labeling
		HS-A5019	E. coli	Coating
		HS-A5194	E. coli	Coating, Labeling
	TP mAb	HS-A5017	Mouse	Labeling



04. Infectious Diseases

Analyte	Product *	Cat. #	Source	Applications in LFA
Tuberculosis (TB)	Recombinant MTB 38KD ESAT6 CFP10 Ag	HS-A5040	E. coli	Coating, Labeling
	Recombinant TB 38 Ag	HS-A5128	E. coli	Coating, Labeling
	Recombinant MTB Mtb8.4 MPT64 TB16.3 Mtb8 Ag	HS-A5129	E. coli	Coating, Labeling
	Recombinant TB MPT64 Ag	HS-A5130	E. coli	Coating, Labeling
	IFN γ mAb	HS-A5126	Mouse	Labeling
		HS-A5127	Mouse	Coating
		HS-A5131	Mouse	Labeling
Typhoid	Recombinant Typhoid Ag	HS-A5020	E. coli	Coating, Labeling



05. Inflammation

Analyte	Product *	Cat. #	Source	Applications in LFA
C-Reaction Protein (CRP)	CRP Ag	HS-E1030	Native	Control Material
	CRP mAb	HS-E1001	Mouse	Coating
		HS-E1002	Mouse	Labeling
		HS-E1011	Mouse	Labeling
		HS-E1012	Mouse	—
		HS-E1013	Mouse	—
		HS-E1014	Mouse	—
		HS-E1015	Mouse	—
		HS-E1016	Mouse	—
		HS-E1018	Mouse	Coating
		HS-E1019	Mouse	Labeling
		HS-E1020	Mouse	Labeling
		HS-E1041	Mouse	—



05. Inflammation

Analyte	Product *	Cat. #	Source	Applications in LFA
Heparin-Binding Protein (HBP)	HBP mAb	HS-E1048	Mouse	Labeling
		HS-E1049	Mouse	Coating
		HS-E1050	Mouse	Labeling
		HS-E1051	Mouse	Coating
Interleukin-6 (IL-6)	Recombinant IL-6 Ag	HS-E1036	E. coli	Control Material
	IL-6 mAb	HS-E1009	Mouse	Labeling
		HS-E1010	Mouse	Coating
		HS-E1024	Mouse	Labeling
		HS-E1025	Mouse	Coating
		HS-E1045	Mouse	Labeling
		HS-E1046	Mouse	Coating
Myxovirusresistance-protein A (MXA)	MXA mAb	HS-E1042	Mouse	—
		HS-E1043	Mouse	—
Procalcitonin (PCT)	Recombinant PCT Ag	HS-E1023	E. coli	Control Material
	PCT mAb	HS-E1007	Mouse	Labeling
		HS-E1008	Mouse	Coating
		HS-E1021	Mouse	Labeling
		HS-E1022	Mouse	Coating
		HS-E1054	Mouse	Labeling
		HS-E1055	Mouse	Coating
		HS-E1056	Mouse	Labeling
HS-E1057	Mouse	Coating		
Serum Amyloid A (SAA)	SAA Ag	HS-E1039	Native	Control Material
	Recombinant SAA Ag	HS-E1017	E. coli	Control Material
	Recombinant SAA Ag	HS-E1026	E. coli	Control Material
	SAA mAb	HS-E1005	Mouse	Labeling
		HS-E1006	Mouse	Coating
		HS-E1027	Mouse	Labeling
		HS-E1028	Mouse	Coating
HS-E1034		Mouse	—	
HS-E1035	Mouse	—		



06. Glucose Metabolism

Analyte	Product *	Cat. #	Source	Applications in LFA
Glutamic Acid Decarboxylase (GAD65)	Recombinant GAD65 Ag	HS-M1014	E. coli	Labeling
		HS-M1032	E. coli	Labeling
Hemoglobin A1C (HbA1c)	HbA1c mAb	HS-G1001	Mouse	—
		HS-G1007	Mouse	Labeling
		HS-G1010	Mouse	Labeling
Hemoglobin (HGB)	HGB mAb	HS-G1002	Mouse	Coating
Protein Tyrosine Phosphatase Antibody (IA-2)	Recombinant IA-2 Ag	HS-M1015	E. coli	Labeling
Zinc Transporter 8 (ZnT8)	Recombinant ZnT8 Ag	HS-M1016	E. coli	Labeling
		HS-M1030	E. coli	Labeling



07. Endocrine

Analyte	Product *	Cat. #	Source	Applications in LFA
Anti-Mullerian Hormone (AMH)	AMH mAb	HS-A1015	Mouse	Labeling
		HS-A1016	Mouse	Coating
Human Chorionic Gonadotropin (HCG)	HCG mAb	HS-A1001	Mouse	Labeling
		HS-A1002	Mouse	Coating
Luteinizing Hormone (LH)	LH mAb	HS-A1003	Mouse	Labeling
		HS-A1004	Mouse	Coating
Progesterone (PROG)	PROG mAb	HS-A1012	Mouse	Labeling
	PROG-BSA Ag	HS-A1013	Synthetic	Coating
Thyroglobulin (TG)	TG mAb	HS-A1028	Mouse	Labeling
		HS-A1029	Mouse	Coating



08. Human Serum

Analyte	Product *	Cat. #	Source	Applications in LFA
Mouse Anti-Human IgG	Mouse anti-human IgG, mAb	HS-A6002	Mouse	Coating, Labeling
		HS-A6003	Mouse	Coating, Labeling
Mouse Anti-Human IgM	Mouse anti-human IgM, mAb	HS-A6005	Mouse	Coating, Labeling
		HS-A6006	Mouse	Coating, Labeling
Anti-Human Serum Albumin (hALB)	hALB pAb	hALBPAB-HS-001	Goat	Labeling
Mouse Anti-Human IgE	Mouse anti-human IgE, mAb	HS-D2002	Mouse	Labeling
		HS-D2003	Mouse	Coating
Anti-Nucleosome	Recombinant Nucleosomes Ag	HS-D2001	E. coli	—



09. Secondary Antibodies

Analyte	Product *	Cat. #	Source	Applications in LFA
Chicken IgY	Chicken IgY, pAb	HS-A6010	Chicken	Labeling
Human IgG	Human IgG, pAb	HS-A6028	Human	—
Human IgM	Human IgM, pAb	HS-A6029	Human	—
Mouse IgG	Mouse IgG, pAb	HS-A6025	Mouse	—
Rabbit IgG	Rabbit IgG, pAb	HS-A6011	Rabbit	Labeling
Goat Anti-Chicken IgY	Goat anti-Chicken IgY, pAb	HS-A6009	Goat	Coating
Goat Anti-Mouse IgG	Goat anti-Mouse IgG, pAb	HS-A6007	Goat	Coating
Goat Anti-Rabbit IgG	Goat anti-Rabbit IgG, pAb	HS-A6008	Goat	Coating
Recombinant Protein A/G	Recombinant Protein A/G Ag	PROAG-HS-001	E. coli	Labeling



10. Animal Diseases

Analyte	Product *	Cat. #	Source	Applications in LFA
Canine Coronavirus (CCV)	CCV mAb	HS-A7001	Mouse	Labeling
		HS-A7002	Mouse	Coating
Canine C-reactive Protein (cCRP)	cCRP Ag	HS-A7015	Native	Control Material
	cCRP mAb	HS-A7018	Mouse	Labeling
		HS-A7019	Mouse	Coating
		HS-A7020	Mouse	Coating
Canine Distemper Virus (CDV)	CDV mAb	HS-A7003	Mouse	Labeling
		HS-A7004	Mouse	Coating
Canine Influenza Virus (CIV)	CIV mAb	HS-A7055	Mouse	Labeling
		HS-A7056	Mouse	Coating
Canine Pancreatic Lipase (CPL)	CPL mAb	HS-A7057	Mouse	Labeling
		HS-A7058	Mouse	Coating
Canine Parvovirus (CPV)	CPV mAb	HS-A7005	Mouse	Labeling
		HS-A7006	Mouse	Coating
Canine Progesterone (cProg)	cProg mAb	HS-A7022	Mouse	Labeling
	cProg-BSA Ag	HS-A7023	Synthetic	Coating
Cat N-Terminal Forebrain Natriuretic Peptide (fNT-proBNP)	fNT-proBNP mAb	HS-A7059	Mouse	Labeling
		HS-A7060	Mouse	Coating
Feline Calicivirus (FCV)	FCV mAb	HS-A7061	Mouse	Labeling
		HS-A7062	Mouse	Coating
Feline Coronavirus (FCoV)	FCOV mAb	HS-A7016	Mouse	Labeling
		HS-A7017	Mouse	Coating
Feline Herpes virus (FHV)	FHV mAb	HS-A7036	Mouse	Labeling
		HS-A7037	Mouse	Coating
Feline Parvovirus (FPV)	FPV mAb	HS-A7007	Mouse	Labeling
		HS-A7008	Mouse	Coating
Feline Serum Amyloid A (fSAA)	fSAA Ag	HS-A7042	Native	Control Material
	fSAA mAb	HS-A7011	Mouse	Labeling
		HS-A7012	Mouse	Coating



10. Animal Diseases

Analyte	Product *	Cat. #	Source	Applications in LFA
Foot-and-Mouth Disease (FMD)	FMD A mAb	HS-A7051	Mouse	Labeling
		HS-A7052	Mouse	Coating
	FMD O mAb	HS-A7053	Mouse	Labeling
		HS-A7054	Mouse	Coating
Rabies Virus (RV)	Recombinant RV-G Ag	HS-M1012	E. coli	Coating, Labeling
	Recombinant RV-N Ag	HS-M1013	E. coli	Coating, Labeling
Secondary Antibody	Cat IgG, pAb	HS-A7045	Cat	Control Material
	Cat IgM, pAb	HS-A7046	Cat	Control Material
	Dog IgG, pAb	HS-A7043	Dog	Control Material
	Dog IgM, pAb	HS-A7044	Dog	Control Material
	Rabbit anti-Cat IgG, pAb	HS-A7009	Rabbit	Coating
	Rabbit anti-Dog IgG, pAb	HS-A7010	Rabbit	Coating
Toxoplasma (TOXO)	Recombinant TOXO Ag	HS-A7031	E. coli	Coating, Labeling

ANIMAL DISEASES





11. Blocking Reagents

Analyte	Product *	Cat. #	Source	Applications in LFA
Antagonist	Antagonist, mAb	HS-A6015	Mouse	—
Heterophilic Antibody	Heterophilic mAb	HS-A6016	Mouse	—
		HS-A6057	Mouse	—
IgG Antagonist	IgG antagonist, mAb	HS-A6033	Mouse	—
IgM Antagonist	IgM antagonist, mAb	HS-A6036	Mouse	—
RBC Antibody	RBC mAb	HS-A6013	Mouse	—
		HS-A6014	Mouse	—
		HS-A6058	Mouse	—



12. Miscellaneous

Analyte	Product *	Cat. #	Source	Applications in LFA
Human P-tau-181 mAb	Human P-tau-181 mAb	HS-A8003	Mouse	Labeling
		HS-A8004	Rabbit	Coating
Procollagen III N-Terminal Propeptide	Procollagen IIINP mAb	HS-G1009	Mouse	Labeling
		HS-G1008	Mouse	Coating
Avidin Reagent	Streptavidin conjugate	HS-D1013	—	Coating
		HS-D1087	—	Coating
	Recombinant Streptavidin Ag	HS-D1003	E. coli	Labeling
Human Cytomegalovirus (CMV)	Recombinant CMV Ag	HS-A1023	E. coli	Coating
Human Single-Cell Virus I (HSV I)	Recombinant HSV I Ag	HS-A1024	E. coli	Coating
Human Single-Cell Virus II (HSV II)	Recombinant HSV II Ag	HS-A1025	E. coli	Coating
Rubella Virus (RV)	Recombinant RV Ag	HS-A1026	E. coli	Coating
Sperm (SP10)	Recombinant SP10 Ag	HS-A1014	E. coli	Control Material
Toxoplasma (TOXO)	Recombinant TOXO Ag	HS-A1011	E. coli	Coating



Healstone Biotech

Healstone Biotech Inc.

Unit #650-655 West Kent Avenue North

Vancouver, B.C., Canada, V6P 6T7

Tel: 1-800-919-3622

1-604-800-2434

info.healstone@healgen.com

www.healstone.com

