

Universal Pipette Tips

- * Compatible with most pipette brands
- * Sterile, RNase, DNase and pyrogen free
- * Autoclavable
- * Low binding option available



HSPT-10P



HSPT-200P



HSPT-1000P



HSPT-10B



HSPT-200B




HSPT-1000B

Cat. No.	P/N	Material	Color	Packaging	Sterile
HSPT-10P	Pipette Tip 10 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	E.O. or E-Beam
HSPT-200P	Pipette Tip 200 μ L	PP	Transparent/Yellow	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSPT-1000P	Pipette Tip 1000 μ L	PP	Transparent/Blue	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSPT-10B	Pipette Tip 10 μ L	PP	Transparent	96 tips/insert, 10 inserts/pk, 10 pks/cs	
HSPT-200B	Pipette Tip 200 μ L	PP	Transparent/Yellow	96 tips/insert, 10 inserts/pk, 10 pks/cs	
HSPT-1000B	Pipette Tip 1000 μ L	PP	Transparent/Blue	96 tips/insert, 5 inserts/pk, 10 pks/cs	

Filter Pipette Tips

- * The filter element prevents contamination by blocking liquids
- * Extended tip
- * Low binding reduces reagent loss
- * Sterile, RNase, DNase and pyrogen free



Cat. No.	P/N	Material	Color	Packaging	Sterile
HSPT-10	Filter Pipette Tip 10 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	E.O. or E-Beam
HSPT-50	Filter Pipette Tip 100 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSPT-200	Filter Pipette Tip 200 μ L	PP	Transparent/Yellow	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSPT-1000	Filter Pipette Tip 1000 μ L	PP	Transparent/Blue	96 tips/rack, 10 racks/pk, 5 pks/cs	

Automation Tips For Hamilton Workstation


Application areas: Genomics, proteomics, cytomics, immunoassays, metabolomics, general liquid handling

* Perfect match with Hamilton Workstations

* Low binding reduces reagent loss

* Transparent/black, filter option available


* Sterile, RNase, DNase and pyrogen free



Cat. No.	P/N	Material	Color	Packaging	Sterile
HSAT-H50	Automation Tip 50 μ L	PP	Black	96 tips/rack, 10 racks/pk, 5 pks/cs	E.O. or E-Beam
HSAT-H300	Automation Tip 300 μ L	PP	Black	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSAT-H1000	Automation Tip 1000 μ L	PP	Black	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSAT-H50A	Automation Tip 50 μ L	PP	Black	96 tips/insert, 4 inserts/pk, 20 pks/cs	
HSAT-H300A	Automation Tip 300 μ L	PP	Black	96 tips/insert, 4 inserts/pk, 20 pks/cs	
HSAT-H1000A	Automation Tip 1000 μ L	PP	Black	96 tips/insert, 4 inserts/pk, 20 pks/cs	
HSAT-H50T	Automation Tip 50 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSAT-H300T	Automation Tip 300 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSAT-H1000T	Automation Tip 1000 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	

Automation Tips For Tecan Workstation


- * Perfect match with Tecan Workstations
- * Low binding reduces reagent loss
- * Transparent/black and filter options available
- * Sterile, RNase, DNase and pyrogen free



Cat. No.	P/N	Material	Color	Packaging	Sterile
HSAT-T50	Automation Tip 50 μ L	PP	Black	96 tips/rack, 10 racks/pk, 5 pks/cs	E.O. or E-Beam
HSAT-T200	Automation Tip 200 μ L	PP	Black	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSAT-T1000	Automation Tip 1000 μ L	PP	Black	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSAT-T50T	Automation Tip 50 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSAT-T200T	Automation Tip 200 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	
HSAT-T1000T	Automation Tip 1000 μ L	PP	Transparent	96 tips/rack, 10 racks/pk, 5 pks/cs	

5 mL & 10 mL Pipette Tips


- * Compatible with Eppendorf, Thermo Fisher pipettes and more (see below for details)
- * Filter option available
- * Pre-filled racks and bulk bags available
- * Sterile, RNase, DNase and pyrogen free



Cat. No.	P/N	Material	Color	Packaging	Sterile
HSPT-5T1	5 mL Pipette tips, compatible with Thermo Finnpiquette, Brand	PP	Transparent	24 tips/rack, 30 racks/cs	E.O. or E-Beam
HSPT-5T2	5 mL Pipette tips, compatible with Thermo Finnpiquette, Brand	PP	Transparent	100 tips/bag, 20 bags/cs	
HSPT-5T3	5 mL Pipette tips, compatible with Eppendorf, Sartorius, DLAB	PP	Transparent	24 tips/rack, 30 racks/cs	
HSPT-5T4	5 mL Pipette tips, compatible with Eppendorf, Sartorius, DLAB	PP	Transparent	100 tips/bag, 20 bags/cs	
HSPT-10T1	10 mL Pipette tips, compatible with Thermo, Finnpiquette	PP	Transparent	24 tips/rack, 30 racks/cs	
HSPT-10T2	10 mL Pipette tips, compatible with Thermo, Finnpiquette	PP	Transparent	100 tips/bag, 20 bags/cs	
HSPT-10T3	10 mL Pipette tips, compatible with Eppendorf, Gilson, Brand, DLAB	PP	Transparent	24 tips/rack, 30 racks/cs	
HSPT-10T4	10 mL Pipette tips, compatible with Eppendorf, Gilson, Brand, DLAB	PP	Transparent	100 tips/bag, 20 bags/cs	

Serological Pipettes


- * Made of high-quality GPPS plastic
- * Individually packed
- * Clear and accurate graduated markings
- * Smooth inner wall for low binding effect
- * Sterile, RNase, DNase and pyrogen free option available



Cat. No.	P/N	Material	Color	Packaging	Sterile
HSSP-1	1 mL Serological Pipette	GPPS	Transparent	50 pcs/bag, 20 bags/cs	E.O. or E-Beam
HSSP-2	2 mL Serological Pipette	GPPS	Transparent	50 pcs/bag, 20 bags/cs	
HSSP-5	5 mL Serological Pipette	GPPS	Transparent	50 pcs/bag, 20 bags/cs	
HSSP-10	10 mL Serological Pipette	GPPS	Transparent	50 pcs/bag, 20 bags/cs	
HSSP-25	25 mL Serological Pipette	GPPS	Transparent	25 pcs/bag, 20 bags/cs	
HSSP-50	50 mL Serological Pipette	GPPS	Transparent	25 pcs/bag, 20 bags/cs	
HSSP-100	100 mL Serological Pipette	GPPS	Transparent	25 pcs/bag, 20 bags/cs	

Pasteur Pipettes

- * Made of transparent polymer polyethylene (PE)
- * Flexible, suitable for micro and uniquely shaped containers
- * With graduated markings
- * Individually packed
- * Sterile, RNase, DNase and pyrogen free



Cat. No.	P/N	Material	Color	Packaging	Sterile
HSPD-1	1 mL Pasteur Pipettes	LDPE	Transparent	500 pcs/bag, 10 bags/cs	E.O. or E-Beam
HSPD-2	2 mL Pasteur Pipettes	LDPE	Transparent	500 pcs/bag, 10 bags/cs	
HSPD-3	3 mL Pasteur Pipettes	LDPE	Transparent	500 pcs/bag, 10 bags/cs	
HSPD-5	5 mL Pasteur Pipettes	LDPE	Transparent	500 pcs/bag, 10 bags/cs	
HSPD-10	10 mL Pasteur Pipettes	LDPE	Transparent	200 pcs/bag, 10 bags/cs	