


PCR Well Plates

Used in molecular cloning, gene expression analysis, genotyping, sequencing and mutagenesis

- * Compatible with most PCR machines
- * High-contrast, black alphanumeric indexing for easy well identification
- * Crystal-clear wells to easily verify sample volume
- * Ultra-thin, uniform wells ensure optimal heat transfer for high reaction efficiency
- * Well ridges reduce sample evaporation when sealing with film, foil, or strip caps
- * Non-skirted, semi-skirted and full-skirted plates with barcode available
- * Sterile, RNase, DNase and pyrogen free



Cat. No.	P/N	Material	Color	Packaging	Sterile
HSPCR-P1	0.1 mL PCR Plate Non-skirted	PP	Transparent	10 pcs/bag, 5 bags/cs	E.O.
HSPCR-P2	0.1 mL PCR Plate Semi-Skirted	PP	Transparent		
HSPCR-P3	0.1 mL PCR Plate Full-Skirted	PP	Transparent		
HSPCR-P4	0.1 mL PCR Plate Non-skirted	PP	White		
HSPCR-P5	0.1 mL PCR Plate Semi-Skirted	PP	White		
HSPCR-P6	0.1 mL PCR Plate Full-Skirted	PP	White		
HSPCR-P7	0.2 mL PCR Plate Non-skirted	PP	Transparent		E-Beam
HSPCR-P8	0.2 mL PCR Plate Semi-Skirted	PP	Transparent		
HSPCR-P9	0.2 mL PCR Plate Full-Skirted	PP	Transparent		
HSPCR-P10	0.2 mL PCR Plate Non-skirted	PP	White		
HSPCR-P11	0.2 mL PCR Plate Semi-Skirted	PP	White		
HSPCR-P12	0.2 mL PCR Plate Full-Skirted	PP	White		

PCR Tubes & 8 Strip Tubes

- * Ultra-thin tube wall for highly efficient heat conduction
- * Compatible with most PCR machines
- * Able to withstand high temperatures
- * Perfect match with lid for secure sealing
- * White and transparent available
- * Sterile, RNase, DNase and pyrogen free



HSPCR-S1T



HSPCR-T2



HSPCR-S2T



HSPCR-T5



HSPCR-S1W

Cat. No.	Description	Material	Color	Packaging	Sterile
HSPCR-T2	0.2 mL PCR Tube with lid	PP	Transparent	100 pcs/bag, 10 bags/cs	E.O. or E-Beam
HSPCR-T5	0.5 mL PCR Tube with lid	PP	Transparent	100 pcs/bag, 10 bags/cs	
HSPCR-S1T	0.1 mL PCR 8 Strip Tube	PP	Transparent	125 pcs/pack, 10 packs/cs	
HSPCR-S2T	0.2 mL PCR 8 Strip Tube	PP	Transparent	125 pcs/pack, 10 packs/cs	
HSPCR-S1W	0.1 mL PCR 8 Strip Tube	PP	White	125 pcs/pack, 10 packs/cs	
HSPCR-S2W	0.2 mL PCR 8 Strip Tube	PP	White	125 pcs/pack, 10 packs/cs	

Deep Well Plates & Magnetic Rod Sets


- * Suitable for multi-channel pipettes and automated pipetting stations
- * Mirror polished inner well lining, achieving less sample residue
- * Good stability and chemical resistance
- * Excellent flatness to ensure the sealing performance of the heat sealing film
- * Stackable, numeric and letter markings for easy identification
- * Autoclavable
- * Can be frozen to -80°C for storage
- * Sterile, RNase, DNase and pyrogen free



HSPCR-D1



HSPCR-M96



HSPCR-M8

Cat. No.	Description	Material	Specification	Packaging	Sterile
HSPCR-D1	Deep Well Plate	PP	1.6 mL	5 pcs/bag, 10 bags/cs	E.O. or E-Beam
HSPCR-D2	Deep Well Plate	PP	2.2 mL	5 pcs/bag, 10 bags/cs	
HSPCR-M8	Magnetic rod set	PP	8/row	2 pcs/bag, 300 bags/cs	
HSPCR-M96	Magnetic rod set	PP	96/row	15 pcs/bag, 10 bags/cs	