




## 02. Cardiovascular

Analyte	Product *	Cat. #	Source	Applications in LFA
Cardiac Troponin I (cTnI)	cTnI mAb	HS-A3001	Mouse	Coating
		HS-A3002	Mouse	Labeling
		HS-A3024	Mouse	Labeling
		HS-A3025	Mouse	Coating
	hs-cTnI mAb	HS-A3049	Mouse	Labeling
		HS-A3050	Mouse	Coating
		HS-A3051	Mouse	Labeling
Cardiac Troponin T (cTnT)	cTnT mAb	HS-A3016	Mouse	Coating
		HS-A3017	Mouse	Labeling
		HS-A3026	Mouse	Labeling
		HS-A3027	Mouse	Coating
Copeptin (CPP)	Recombinant CPP Ag	HS-M1026	E. coli	Control Material
	CPP mAb	HS-M1020	Mouse	Labeling
		HS-M1021	Mouse	Coating
Creatine Kinase Isoenzyme-MB (CK-MB)	CK-MB mAb	HS-A3007	Mouse	Coating
		HS-A3008	Mouse	Labeling
		HS-A3028	Mouse	Labeling
		HS-A3029	Mouse	Coating
		HS-A3053	Mouse	—
		HS-A3054	Mouse	—
D-Dimer (D-D)	D-D mAb	HS-A3020	Mouse	Coating
		HS-A3021	Mouse	Labeling
		HS-A3032	Mouse	Labeling
		HS-A3033	Mouse	Coating
		HS-A3037	Mouse	Labeling
		HS-A3038	Mouse	Coating





## 02. Cardiovascular

Analyte	Product *	Cat. #	Source	Applications in LFA
Growth Stimulation Expressed Gene 2 (ST2)	Recombinant ST2 Ag	HS-A3046	293 cells	Control Material
	ST2 mAb	HS-A3047	Mouse	Labeling
		HS-A3048	Mouse	Coating
Myoglobin (Myo)	Myo mAb	HS-A3004	Mouse	Coating
		HS-A3005	Mouse	Labeling
		HS-A3030	Mouse	Labeling
		HS-A3031	Mouse	Coating, Labeling
		HS-A3040	Mouse	Coating
N-terminal Brain Natriuretic Peptide Precursor (NT proBNP)	NT proBNP mAb	HS-A3022	Mouse	Coating
		HS-A3023	Mouse	Labeling
		HS-A3034	Mouse	Labeling
		HS-A3035	Mouse	Coating
		HS-A3039	Mouse	Labeling
Oxidized Low Density Lipoprotein (OxLDL)	OxLDL Ag	HS-M1027	Native	Control Material
	OxLDL mAb	HS-M1022	Mouse	Labeling
		HS-M1023	Mouse	Coating
S100 $\beta$ Protein (S100 $\beta$ )	Recombinant S100 $\beta$ Ag	HS-M1028	E. coli	Control Material
		HS-M1036	E. coli	Control Material
	S100 $\beta$ mAb	HS-M1024	Mouse	Labeling
		HS-M1025	Mouse	Coating



Fertility & Women's Health



Drugs



Infectious Diseases



Oncology



Cardiovascular



Autoimmunity



Metabolism