



06. Glucose Metabolism

Analyte	Product *	Cat. #	Source	Applications in LFA
Glutamic Acid Decarboxylase (GAD65)	Recombinant GAD65 Ag	HS-M1014	E. coli	Labeling
		HS-M1032	E. coli	Labeling
Hemoglobin A1C (HbA1c)	HbA1c mAb	HS-G1001	Mouse	—
		HS-G1007	Mouse	Labeling
		HS-G1010	Mouse	Labeling
Hemoglobin (HGB)	HGB mAb	HS-G1002	Mouse	Coating
Protein Tyrosine Phosphatase Antibody (IA-2)	Recombinant IA-2 Ag	HS-M1015	E. coli	Labeling
Zinc Transporter 8 (ZnT8)	Recombinant ZnT8 Ag	HS-M1016	E. coli	Labeling
		HS-M1030	E. coli	Labeling



07. Endocrine

Analyte	Product *	Cat. #	Source	Applications in LFA
Anti-Mullerian Hormone (AMH)	AMH mAb	HS-A1015	Mouse	Labeling
		HS-A1016	Mouse	Coating
Human Chorionic Gonadotropin (HCG)	HCG mAb	HS-A1001	Mouse	Labeling
		HS-A1002	Mouse	Coating
Luteinizing Hormone (LH)	LH mAb	HS-A1003	Mouse	Labeling
		HS-A1004	Mouse	Coating
Progesterone (PROG)	PROG mAb	HS-A1012	Mouse	Labeling
	PROG-BSA Ag	HS-A1013	Synthetic	Coating
Thyroglobulin (TG)	TG mAb	HS-A1028	Mouse	Labeling
		HS-A1029	Mouse	Coating