



## 03. Tumor Markers

Analyte	Product *	Cat. #	Source	Applications in LFA
Alpha-Fetoprotein (AFP)	AFP mAb	HS-A2007	Mouse	Labeling
		HS-A2008	Mouse	Coating
		HS-A2044	Mouse	Labeling
		HS-A2045	Mouse	Coating
		HS-A2054	Mouse	Coating
Carbohydrate Antigen 125 (CA125)	CA125 mAb	HS-A2022	Mouse	Labeling
		HS-A2023	Mouse	Coating
		HS-A2073	Mouse	Labeling
		HS-A2074	Mouse	Coating
Carbohydrate Antigen 153 (CA153)	CA153 mAb	HS-A2047	Mouse	Labeling
		HS-A2048	Mouse	Coating
Carbohydrate Antigen 199 (CA199)	CA199 mAb	HS-A2034	Mouse	Coating
		HS-A2035	Mouse	Labeling
Carbohydrate Antigen 242 (CA242)	CA242 mAb	HS-A2086	Mouse	Labeling
		HS-A2087	Mouse	Coating
Carcinoembryonic Antigen (CEA)	CEA mAb	HS-A2005	Mouse	Labeling
		HS-A2006	Mouse	Coating
Cytokeratin (Cy 211)	Cy 211 mAb	HS-A2029	Mouse	Coating
		HS-A2046	Mouse	Labeling
		HS-A2071	Mouse	Labeling
Fecal Occult Blood (FOB)	FOB mAb	HS-A2001	Mouse	Labeling
		HS-A2002	Mouse	Coating
		HS-A2025	Mouse	Labeling
		HS-A2026	Mouse	Coating
Ferritin (FER)	FER mAb	HS-A2040	Mouse	Labeling
		HS-A2041	Mouse	Coating
Gastrin 17 (G17)	Recombinant G17 Ag	HS-M1029	E. coli	Control Material
	G17 mAb	HS-M1018	Mouse	Labeling
		HS-M1019	Mouse	Coating
		HS-M1034	Mouse	Coating
		HS-M1035	Mouse	Labeling
Human Epididymis Protein 4 (HE4)	HE4 mAb	HS-A2036	Mouse	Labeling
		HS-A2037	Mouse	Coating
		HS-A2060	Mouse	Labeling
		HS-A2085	Mouse	Labeling



## 03. Tumor Markers

Analyte	Product *	Cat. #	Source	Applications in LFA
Neuron-Specific Enolase (NSE)	NSE mAb	HS-A2038	Mouse	Labeling
		HS-A2072	Mouse	Coating
Pepsinogen (PG)	PG I mAb	HS-A2049	Mouse	Labeling
		HS-A2050	Mouse	Coating
		HS-A2061	Mouse	Labeling
		HS-A2062	Mouse	Coating
		HS-A2063	Mouse	Labeling
		HS-A2064	Mouse	Coating
	PG II mAb	HS-A2051	Mouse	Labeling
		HS-A2052	Mouse	Coating
		HS-A2065	Mouse	Labeling
		HS-A2066	Mouse	Coating
		HS-A2067	Mouse	Labeling
		HS-A2068	Mouse	Coating
Pro-Gastrin-Releasing Peptide (ProGRP)	ProGRP mAb	HS-A2027	Mouse	Coating
		HS-A2028	Mouse	Labeling
Prostate Specific Antigen (PSA)	PSA mAb	HS-A2003	Mouse	Labeling
		HS-A2004	Mouse	Coating
Squamous Cell Carcinoma (SCCA)	SCCA mAb	HS-A2019	Mouse	Labeling
		HS-A2020	Mouse	Coating
		HS-A2069	Mouse	Labeling
		HS-A2070	Mouse	Coating
Vascular Endothelial Growth Factor (VEGF)	Recombinant VEGF Ag	HS-A2059	E. coli	Control Material
	VEGF mAb	HS-A2055	Mouse	Labeling
		HS-A2056	Mouse	Labeling
		HS-A2057	Mouse	Labeling
	VEGF pAb	HS-A2058	Goat	Coating
HS-A2096		Mouse	Coating	