

## Disposable Swabs

Suitable for sampling biological specimens in a variety of scenarios

- \* Polypropylene, polyamide, and polyurethane sampling swab tips are available
- \* Swab rods of different lengths and break points to suit a variety of scenarios
- \* Independent packaging, customized available
- \* Sterile, RNase, DNase and pyrogen free



HSDS-0601



HSDS-0701



HSDS-1001



HSDS-1201



HSDS-1301

HSDS-1801

| Cat. No.  | Description   | Application   | Sterile              |
|-----------|---|---|----------------------|
| HSDS-0601 | L80 mm, polyamide 66 sampling tip                     | Transfer or collect biological samples              | E.O.<br>or<br>E-Beam |
| HSDS-0701 | L80 mm, polyurethane sampling tip                     | Collect biological samples from the nose and throat |                      |
| HSDS-0901 | L150 mm, polyamide 66 sampling tip                    | Collect biological samples from the throat          |                      |
| HSDS-1001 | L150 mm, polyurethane sampling tip                    | Collect biological samples from the throat          |                      |
| HSDS-1101 | L150 mm, polyamide 66 sampling tip                    | Collect biological samples from the throat          |                      |
| HSDS-1201 | L150 mm, 80 mm break point, polyamide 66 sampling tip | Collect biological samples from the throat          |                      |
| HSDS-1301 | L150 mm, 80 mm break point, polyurethane sampling tip | Collect biological samples from the throat          |                      |
| HSDS-1401 | L150 mm, 80 mm break point, polyamide 66 sampling tip | Collect biological samples from the nose and throat |                      |
| HSDS-1501 | L150 mm, 40 mm break point, polyurethane sampling tip | Collect biological samples from the throat          |                      |
| HSDS-1601 | L150 mm, 40 mm break point, polyurethane sampling tip | Collect biological samples from the throat          |                      |
| HSDS-1701 | L150 mm, 40 mm break point, polyamide 66 sampling tip | Collect biological samples from the throat          |                      |
| HSDS-1801 | L150 mm, 80 mm break point, polyamide 66 sampling tip | Collect biological samples from the throat          |                      |

## Transport Swabs with Tube

Used to collect microbial samples from skin, wounds, throat, vagina and more, storing them for short periods of time

- \* Each cotton swab and storage tube set wrapped in individual sterilization pouches
- \* Writing area for easy labeling
- \* Custom labeling and packaging available
- \* Sterile, RNase, DNase and pyrogen free



HSTS-1

## Transport Swabs with Medium

- \* Suitable for sampling from pharyngeal, nose, ear, eye, vagina, urethra, skin, intestinal tract, etc
- \* Ideal for the collection, storage, and transportation of microorganism specimens, protecting their character and activity before analysis
- \* Can be stored for 72 hours after collection
- \* Sterile, RNase, DNase and pyrogen free



HSTSM-1



HSTSM-2



HSTSM-3

| Cat. No. | P/N                      | Material | Tube Size    | Packaging               | Sterile              |
|----------|--------------------------|----------|--------------|-------------------------|----------------------|
| HSTS-1   | Transport Swab with Tube | PP       | Ø13 x 175 mm | 100 pcs/bag, 20 bags/cs | E.O.<br>or<br>E-Beam |
| HSTSM-1  | Amies Agar Gel Swab      | PP       | Ø13 x 175 mm | 100 pcs/bag, 20 bags/cs |                      |
| HSTSM-2  | Stuart Agar Gel Swab     | PP       | Ø13 x 175 mm | 100 pcs/bag, 20 bags/cs |                      |
| HSTSM-3  | Cary-Blair Agar Gel Swab | PP       | Ø13 x 175 mm | 100 pcs/bag, 20 bags/cs |                      |