


Test Tubes

Made of superior PET/PS material with good chemical compatibility; suitable for most polar organic solvents, weak acids, and weak bases

- * Maximum force: 4,500 rpm
- * Multiple sizes and types available
- * Custom labeling available
- * Matching caps available to prevent liquid leakage




Cat. No.	P/N	Material	Size	Packaging	Sterile
HSTT-1	Test Tube	PS/PET	Ø10 x 75 mm	500 pcs/bag, 10 bags/cs	E.O. or E-Beam
HSTT-2	Test Tube	PS/PET	Ø12 x 75 mm	500 pcs/bag, 10 bags/cs	
HSTT-3	Test Tube	PS/PET	Ø10 x 100 mm	500 pcs/bag, 10 bags/cs	
HSTT-4	Test Tube	PS/PET	Ø12 x 100 mm	500 pcs/bag, 10 bags/cs	
HSTT-5	Test Tube	PS/PET	Ø13 x 100 mm	500 pcs/bag, 10 bags/cs	
HSTT-6	Test Tube	PS/PET	Ø15 x 100 mm	300 pcs/bag, 10 bags/cs	
HSTT-7	Test Tube	PS/PET	Ø16 x 100 mm	200 pcs/bag, 10 bags/cs	

Urinary Sediment Tubes

Used for pouring collected urine into a test tube for automated instrument analysis

* High transparency for easy observation of sample volume

* Capless, screw cap, and pierceable cap options available



HSUST-1



HSUST-2



HSUST-3

Cat. No.	P/N	Material	Specification	Packaging	Sterile
HSUST-1	Urinary Sediment Tube	PS	12 mL, with screw lid	200 pcs/bag, 10 bags/cs	E.O.
HSUST-2	Urinary Sediment Tube	PS	12 mL, with pierceable lid	200 pcs/bag, 10 bags/cs	or
HSUST-3	Urinary Sediment Tube	PP	12 mL, with pierceable lid	200 pcs/bag, 10 bags/cs	E-Beam