



# Malaria P.f./Pan Ag Rapid Test




This is a rapid lateral flow chromatographic immunoassay for the simultaneous detection and differentiation of HRP-II (Histidine-rich protein II) specific to *Plasmodium falciparum* and pLDH (Plasmodium lactate dehydrogenase) specific to *Plasmodium* species (Pan) in human blood specimen as an aid in the diagnosis of Malaria infection.

## FEATURES

- Health Canada Class II licensed
- Easy to use
- Read results in 20–30 minutes
- Test can be done with finger prick blood

## SPECIFICATIONS

- Individually packaged cassette: 25 tests/kit
- Storage: 2-30°C/36-86°F
- Shelf life: 24 months



| Product Description                                   | Catalog No.        | Kit Size     |
|-------------------------------------------------------|--------------------|--------------|
| Malaria P.F./Pan Ag Rapid Test Cassette (Whole Blood) | GCMAL(pf/pan)-402a | 25 Tests/Kit |



# RSV Antigen Test



The RSV Antigen Rapid Test cassette is rapid chromatographic immunoassay for the detection of respiratory syncytial virus (RSV) antigen (viral fusion protein) in human nasal/nasopharyngeal swab, nasal/nasopharyngeal aspirate samples. It is intended as a screening test and as an aid in the diagnosis of RSV respiratory infections.

## FEATURES

- Health Canada Class II licensed.
- Rapid and accurate results within 15 minutes.
- Applicable for both nasal/nasopharyngeal swab, and nasal suction fluid use

## SPECIFICATIONS

- 20 tests/kit
- Storage: 2-30°C/36-86°F
- Shelf life: 24 months



| Product Description                    | Catalog No. | Kit Size     |
|----------------------------------------|-------------|--------------|
| RSV Antigen Rapid Test Cassette (Swab) | GCRSV-502a  | 20 Tests/Kit |