



# Harm Reduction



## Fentanyl/Xylazine Residue Test Strip


A rapid screening test for the qualitative detection of Fentanyl or Xylazine residue (e.g. dried liquid or powder) on surfaces, including the surface of pills.

This test is used to prevent misuse and overdose deaths of Fentanyl and Xylazine.

### FEATURES

- All necessary tools (swabs and buffers) are included in the test kits
- Does not require the crushing of pills and/or dissolving powders to liquids before testing
- Easy collecting and testing procedures
- Easy visual interpretation of results
- Wait time to results: 3-10 minutes  
(Do not interpret results after 1 hour)
- Fentanyl/Xylazine can be detected at 5ng/ml and 50 ng/ml respectively. More sensitive in detection windows than other products
- Most FEN derivatives can be detected with low cut off values
- 250+ analytes were tested and shown in both IFU
- This is a non-IVD product

### Test Procedure - Strip



FEN test strip



XYL test strip



Product Description	Catalog No.	Format	Cut-off Value (ng/mL)	Kit Size
Fentanyl Residue Test Strip	GBFEN-1001a	Strip	5	20 Tests/Kit
Xylazine Residue Test Strip	GBXYL-1001a	Strip	50	20 Tests/Kit



# Harm Reduction



## Fentanyl/Xylazine Residue Test Cassette Combo


A rapid screening test for the qualitative detection of Fentanyl or Xylazine residue (e.g. dried liquid or powder) on surfaces, including the surface of pills.

This test is used to prevent misuse and overdose deaths of Fentanyl and Xylazine.

### FEATURES

- All necessary tools (integrated swabs and buffers) are included in the test kits
- Does not require the crushing of pills and/or dissolving powders to liquids before testing
- Easy collecting and testing procedures
- Easy visual interpretation of results
- Wait time to results: 5-10 minutes (Do not interpret results after 1 hour)
- Fentanyl/Xylazine can be detected at 5ng/ml and 50 ng/ml respectively. More sensitive in detection windows than other products
- Most FEN derivatives can be detected with low cut off values
- 250+ analytes were tested and shown in both IFU
- This is a non-IVD product

### Test Procedure - Cassette



Product Description	Catalog No.	Format	Cut-off Value (ng/mL)	Kit Size
Fentanyl/Xylazine Residue Test Cassette	GBDOA-1025-001	Cassette	5 (FEN) 50 (XYL)	20 Tests/Kit

**Disclaimer:** Healgen/Healstone does not encourage the use, supply, or production of illegal drugs or controlled substances in any way. The device is intended for harm reduction purposes. Follow the advice of your local harm reduction or public health agency. A negative result does not at any time rule out the presence of FEN/XYL as they may be present below the minimum detection level of the test.