



Urinary Tract Infection Test Strips



The URcheck™ URS Test Strip is for the qualitative detection of leukocytes (white blood cells) and Nitrite in urine as an aid in the screening of a urinary tract infection (UTI). It is intended for near patient and over-the-counter, home use only.

Lower urinary tract infections may involve the bladder (cystitis) or urethra (urethritis) while upper urinary tract infections affect the kidney (pyelonephritis).


FEATURES

- Health Canada Class III licensed for home use
- Rapid and accurate results within 2 minutes
- High test sensitivity:
 Nitrite (mg/dL): 0.08
 Leukocytes (cells/ μ L): 15
- Convenient to hold and easy to test with the thumb handler
- English and French bilingual labelling



SPECIFICATIONS

- 1/3/6 tests per kit package available
- Storage: 15-30°C/59-86°F



Ordering Information

Product Description	Catalog No.	Kit Size
URcheck™ URS Test Strip	UTI-2T	1/3/6 tests/kit