



Urine Collection Containers



FEATURES

- Made with polypropylene allowing for clarity and durability
- Data supported leak-proof capacity
- Gradations on cup exterior for easy reading
- High versatility: Ideal for strips, cassettes, and 1 to 12 panel drug tests

SPECIFICATIONS

- 90, 120, 130 ml volume options
- Sterile cup equipped with and/or without the temperature strip and tamper stick
- HSUC-5/HSUC-5H: Taper cup, convenient for collection. Audible click indicates secure capping.
- HSUC-7/HSUC-7H: Audible click indicates secure capping. Cap with O-ring ensures secure sealing during transport to laboratories for confirmatory testing



Cat. #	Description	Dimensions	Volume	O Ring	Leak proof pressure	Packaging
HSUC-3	90 mL straight urine container	Ø 47 X 69 mm	90 mL	No	96 Kpa	100 /carton
HSUC-3H	90 mL straight urine container with temperature strip and tamper sticker	Ø 47 X 69 mm	90 mL	No	96 Kpa	100 /carton
HSUC-5	120 mL taper urine container with clicker	Ø 68/50 X 79mm	120 mL	No	200 Kpa	100 /carton
HSUC-5H	120 mL taper urine container with clicker, temperature strip and tamper sticker	Ø 68/50 X 79mm	120 mL	No	200 Kpa	100 /carton
HSUC-7	130 mL straight urine container with clicker	Ø 57 X 83 mm	130 mL	Yes	200 Kpa	100 /carton
HSUC-7H	130 mL straight urine container with clicker, temperature strip and tamper sticker	Ø 57 X 83 mm	130 mL	Yes	200 Kpa	100 /carton