



Harm Reduction




Residue Test Strip


A rapid screening test for the qualitative detection of Fentanyl, Xylazine, Benzodiazepine, Carfentanyl, Medetomidine, Nitazene residue (e.g. dried liquid or powder) on surfaces, including the surface of pills. This test is used to prevent misuse and overdose deaths of extremely harmful and dangerous drugs.

FEATURES


- All necessary tools (swabs and buffers) are included in the test kits
- Does not require the crushing of pills and/or dissolving powders to liquids before testing
- Easy collecting and testing procedures
- Easy visual interpretation of results
- Wait time to results: 3-10 minutes
- More sensitive in detection windows than other products
- Most FEN derivatives can be detected with low cut off values
- Carfentanyl has no cross reaction with Fentanyl at < 50,000 ng/ml
- Fentanyl, Nitazene have no cross reactions with Caffeine
- Up to 3 test items in 1 test strip
- This is a non-IVD product



Test strips



Product Description	Catalog No.	Format	Cut-off Value (ng/mL)	Kit Size
Fentanyl Test Strip	GBFEN-1001a	Strip	5	3, 5, 20 Tests/Kit
Xylazine Test Strip	GBXYL-1001a	Strip	50	3, 5, 20 Tests/Kit
Benzodiazepine Test Strip	GBBZO-1001a	Strip	10	3, 5, 20 Tests/Kit
Carfentanyl Test Strip	TBC	Strip	1000	3, 5, 20 Tests/Kit
Medetomidine Test Strip	GBMDT-1001b	Strip	100	3, 5, 20 Tests/Kit
Nitazene Test Strip	GBNTZ-1001Ba	Strip	100	3, 5, 20 Tests/Kit
FEN5/XYL50 Combo Test Strip	TBC	Strip	5 (FEN), 50 (XYL)	3, 5, 20 Tests/Kit
FEN5/XYL50/BZO10 Combo Test Strip	TBC	Strip	5 (FEN) 50 (XYL), 10 (BZO)	3, 5, 20 Tests/Kit



Harm Reduction




Multi Panel Drugs Test Combo Cassette

A rapid screening test for the qualitative detection of Fentanyl, Xylazine, Benzodiazepine, Carfentanyl, Medetomidine, Nitazene residue (e.g. dried liquid or powder) on surfaces, including the surface of pills. This test is used to prevent misuse and overdose deaths of extremely harmful and dangerous drugs.

FEATURES

- All necessary tools (integrated swabs and buffers) are included in the test kits
- Does not require the crushing of pills and/or dissolving powders to liquids before testing
- Easy collecting and testing procedures
- Easy visual interpretation of results
- Wait time to results: 5-10 minutes
- More sensitive in detection windows than other products
- Most FEN derivatives can be detected with low cut off values
- Carfentanyl has no cross reaction with Fentanyl at < 50,000 ng/ml
- Fentanyl, Nitazene have no cross reactions with Caffeine
- Up to 5 test items in 1 cassette
- This is a non-IVD product



Product Description (Combo Test Cassettes)	Catalog No.	Cut-off Value (ng/mL)	Kit Size
FEN5/XYL50	GBDOA-1025-001	5 (FEN), 50 (XYL)	3, 5, 20 Tests/Kit
FEN5/XYL50/BZO10	TBC	5 (FEN), 50 (XYL) 10 (BZO)	3, 5, 20 Tests/Kit
FEN5/CFN1000	GBDOA-1025-002	5 (FEN), 1000 (CFN)	3, 5, 20 Tests/Kit
FEN5/XYL50/BZO10/MDT100	TBC	5 (FEN), 50 (XYL), 10 (BZO), 100 (MDT)	3, 5, 20 Tests/Kit
FEN5/CFN1000/XYL50/BZO10/MDT100	TBC	5 (FEN), 1000 (CFN), 50 (XYL) 10 (BZO), 100 (MDT)	3, 5, 20 Tests/Kit

Disclaimer: Healgen/Healstone does not encourage the use, supply, or production of illegal drugs or controlled substances in any way. The device is intended for harm reduction purposes. Follow the advice of your local harm reduction or public health agency. A negative result does not at any time rule out the presence of FEN/XYL as they may be present below the minimum detection level of the test.