



Malaria P.f./Pan Ag Rapid Test




This is a rapid lateral flow chromatographic immunoassay for the simultaneous detection and differentiation of HRP-II (Histidine-rich protein II) specific to *Plasmodium falciparum* and pLDH (Plasmodium lactate dehydrogenase) specific to *Plasmodium* species (Pan) in human blood specimen as an aid in the diagnosis of Malaria infection.

FEATURES

- Health Canada Class II licensed
- Easy to use
- Read results in 20–30 minutes
- Test can be done with finger prick blood

SPECIFICATIONS

- Individually packaged cassette: 25 tests/kit
- Storage: 2-30°C/36-86°F
- Shelf life: 24 months



Product Description	Catalog No.	Kit Size
Malaria P.F./Pan Ag Rapid Test Cassette (Whole Blood)	GCMAL(pf/pan)-402a	25 Tests/Kit